

MSI BY PCR ANALYSIS: MP24-190MP00053

Order: 9999999999 - Reflex for Order 0000000000

Collected MM/DD/YYYY 0:00 PM Status: Final result Visible to patient: Yes (seen)

Dx: Person encountering health services t...

Component

Microsatellite Status**PCR Analysis**

Cleveland Clinic MSI PCR Analysis

Laboratory Accession Number: 0000000000

Case #: P00-000000

Part ID: XXXXXX

Sample Type: FFPET

% Tumor: 15

RESULT:

Indeterminate

INTERPRETATION:

This sample is indeterminate for MSI testing by PCR. Microsatellite instability is NOT DETECTED in this specimen; however, the tumor purity was significantly less than required (<40%). As such, a negative result may be due to decreased sensitivity and is considered indeterminate. Testing of an alternative specimen is recommended. Clinical and pathological correlation required.

Microsatellite instability (MSI) is a hypermutator phenotype, caused by the loss of mismatch repair (MMR) gene function (MMR deficiency; dMMR). Tumor MSI-H/dMMR status can be due to germline genetic mutation (Lynch Syndrome) or, more commonly, somatic genetic and epigenetic changes. MSI/dMMR status is routinely assessed for any sporadic cancer type in the Lynch Syndrome spectrum (i.e., colorectal, endometrial, etc.), and MSI-H/dMMR status can predict clinical benefit of immune checkpoint blockade for solid tumors. Specifically, the anti-PD-1 antibody pembrolizumab is an FDA-approved drug for therapy of adult and pediatric patients with unresectable or metastatic MSI-H/dMMR solid cancers that have progressed following prior treatment. In addition, the anti-PD-1 antibody dostarlimab in combination with carboplatin and paclitaxel is FDA-approved for the treatment of patients with MSI-H/dMMR endometrial cancer. The anti-PD-1 antibodies pembrolizumab or nivolumab, as single-agents, and the anti-CTLA4 antibody ipilimumab in combination with nivolumab are FDA-approved for the treatment of patients with MSI-H/dMMR metastatic colorectal cancer.

METHODS:

The Idylla MSI PCR analysis assay is a rapid, automated, cartridge-based system to detect microsatellite instability in seven monomorphic DNA biomarkers (ACVR2A, BTBD7, DIDO1, MRE11, RYR3, SEC31A and SULF2) using real-time PCR from formalin-fixed paraffin embedded tissue specimens. The assay was performed according to manufacturer's instructions. Samples with less than 40 percent tumor purity undergo

macrodissection prior to analysis. This assay is performed on the Idylla system (Biocartis US, Itasca, IL). MSI PCR cartridges contain dried reagents (enzymes, primers, and probes), 1.7 ml liquefaction buffer and 2.2 ml dilution buffer to enable tissue to undergo liquefaction, cell lysis, extraction, real-time PCR amplification and detection. After PCR amplification, molecular beacons differentially dissociate from the wild type or mutated amplicons with increasing temperature. The MSI specific software included in the Idylla system checks the validity of the measured fluorescence profiles. These profiles are analyzed by the software to determine marker genotype. The MSI status of the sample can be determined if at least five valid marker-specific fluorescence profiles could be fully analyzed (otherwise the MSI status will be called "Invalid"). If the sample is valid for test resulting, at least two mutant markers will result in a status being MSI-H (Microsatellite Instability-High), otherwise the status will be scored as MSS (Microsatellite Stable).

LIMITATIONS:

Specimens must contain at least 40 percent tumor cells; if less than 40 percent tumor is utilized, a negative result is of indeterminate significance. Tumor heterogeneity, tumor burden, specimen degradation or other limitations of the technology may affect the sensitivity. Interfering substances, specifically formalin, decalcification agents, fixation agents containing heavy metals or preservation of buffy coats using Hank's Balanced Salt Solution (HBSS) can potentially affect assay performance. This test cannot distinguish between somatic and germline variants. Molecular microsatellite instability can be assessed by clinically validated IHC, PCR and NGS methods; discrepancies between these results can occur due to technical or biological reasons. Test results obtained using this assay must be interpreted by healthcare professionals in conjunction with other clinical findings, family history and other laboratory data.

REFERENCES:

- 1) Gupta R et al. The impact of microsatellite stability status in colorectal cancer. *Curr Probl Cancer*. 2018 Nov;42(6):548-559. doi: 10.1016/j.currproblcancer.2018.06.010. Epub 2018 Jul 18. PMID: 30119911.
- 2) Luchini C et al. ESMO recommendations on microsatellite instability testing for immunotherapy in cancer, and its relationship with PD-1/PD-L1 expression and tumour mutational burden: a systematic review-based approach. *Ann Oncol*. 2019 Aug 1;30(8):1232-1243. doi: 10.1093/annonc/mdz116. PMID: 31056702.
- 3) Hause RJ et al. Classification and characterization of microsatellite instability across 18 cancer types. *Nat Med*. 2016 Nov;22(11):1342-1350. doi: 10.1038/nm.4191. Erratum in: *Nat Med*. 2017 Oct 6;23(10):1241. Erratum in: *Nat Med*. 2018 Apr 10;24(4):525. PMID: 27694933.
- 4) Zwaenepoel K et al. Clinical Performance of the Idylla MSI Test for a Rapid Assessment of the DNA Microsatellite Status in Human Colorectal Cancer. *J Mol Diagn*. 2020 Mar;22(3):386-395. doi: 10.1016/j.moldx.2019.12.002. Epub 2019 Dec 24. PMID: 31881332.
- 5) Ukkola I et al. Detection of microsatellite instability with Idylla MSI assay in colorectal and endometrial cancer. *Virchows Arch*. 2021

Sep;479(3):471-479. doi: 10.1007/s00428-021-03082-w. PMID: 33755781; PMCID: PMC8448708.

6) Dedeurwaerdere F et al. Comparison of microsatellite instability detection by immunohistochemistry and molecular techniques in colorectal and endometrial cancer. *Sci Rep.* 2021 Jun 18;11(1):12880. doi: 10.1038/s41598-021-91974-x. PMID: 34145315; PMCID: PMC8448708.

7) Nadorvari ML et al. Comparison of standard mismatch repair deficiency and microsatellite instability tests in a large cancer series. *J Transl Med.* 2024 Feb 13;22(1):150. doi: 10.1186/s12967-024-04960-y. PMID: 38350968; PMCID: PMC10863158.

8) Le DT et al. Mismatch repair deficiency predicts response of solid tumors to PD-1 blockade. *Science.* 2017 Jul 28;357(6349):409-413. doi: 10.1126/science.aan6733. PMID: 28596308; PMCID: PMC5576142.

9) Overman MJ et al. Durable Clinical Benefit With Nivolumab Plus Ipilimumab in DNA Mismatch Repair-Deficient/Microsatellite Instability-High Metastatic Colorectal Cancer. *J Clin Oncol.* 2018 Mar 10;36(8):773-779. doi: 10.1200/JCO.2017.76.9901. PMID: 29355075.

10) Abida W et al. Analysis of the Prevalence of Microsatellite Instability in Prostate Cancer and Response to Immune Checkpoint Blockade. *JAMA Oncol.* 2019 Apr 1;5(4):471-478. doi: 10.1001/jamaoncol.2018.5801. PMID: 30589920; PMCID: PMC6459218.

11) Le DT et al. Phase II Open-Label Study of Pembrolizumab in Treatment-Refractory, Microsatellite Instability-High/Mismatch Repair-Deficient Metastatic Colorectal Cancer: KEYNOTE-164. *J Clin Oncol.* 2020 Jan 1;38(1):11-19. doi: 10.1200/JCO.19.02107. PMID: 31725351; PMCID: PMC7031958.

12) Yoshino T et al. JSCO-ESMO-ASCO-JSMO-TOS: international expert consensus recommendations for tumour-agnostic treatments in patients with solid tumours with microsatellite instability or NTRK fusions. *Ann Oncol.* 2020 Jul;31(7):861-872. doi: 10.1016/j.annonc.2020.03.299.. PMID: 32272210.

13) Suehnholz SP et al. Quantifying the Expanding Landscape of Clinical Actionability for Patients with Cancer. *Cancer Discov.* 2024 Jan 12;14(1): 49-65. doi: 10.1158/2159-8290.CD-23-0467. PMID: 37849038; PMCID: PMC10784742.

14) Mirza MR et al. Dostarlimab for Primary Advanced or Recurrent Endometrial Cancer. *N Engl J Med.* 2023 Jun 8;388(23):2145-2158. doi: 10.1056/NEJMoa2216334. PMID: 36972026.

15) Idylla Operator Manual. Biocartis BCT018459. May2023.

16) Assay Instructions for Use Idylla MSI Test. Biocartis BCT012834 Version 1. February 2023.

DISCLAIMER:

This test was developed and its performance characteristics determined by Cleveland Clinic's Robert J. Tomsich Pathology and Laboratory Medicine Institute (RT-PLMI). It has not been cleared or approved by the FDA. RT-PLMI is regulated under CLIA as certified to perform high-complexity testing. This test is used for clinical purposes. It should not be regarded as investigational or for research.

Testing and interpretation performed at Cleveland Clinic, 9500 Euclid Ave, Cleveland, OH 44195. CLIA Number: 36D0656094

As reviewed by Elizabeth Azzato, MD, PhD