

# High Sensitivity Flow Cytometry for Paroxysmal Nocturnal Hemoglobinuria

## Background Information

Paroxysmal nocturnal hemoglobinuria (PNH) is a clonal stem cell disorder characterized by mutations in the PIGA gene leading to loss of cell surface proteins linked to glycosphosphatidylinositol (GPI) anchors. Patients affected by PNH display complement mediated hemolysis, thrombosis and bone marrow failure though the clinical presentation is variable. The presence of a PNH clone occurs in classical hemolytic PNH, generally at levels above 1%; however, PNH clones may also be seen in other disorders such as aplastic anemia and myelodysplastic syndrome. Flow cytometric immunophenotypic analysis is the method of choice to detect populations of GPI anchor-deficient cells used in the diagnosis of PNH and to monitor patients with an established diagnosis.<sup>1</sup> PNH erythrocyte clones may be divided into those with partial loss of CD59 (PNH Type II Cells) or complete loss of CD59 (PNH Type III Cells).

Cleveland Clinic Laboratories now offers high sensitivity flow cytometry testing for PNH using whole blood. This procedure labels red blood cells (RBC) and white blood cells (WBC) for the detection of GPI-linked surface antigens using monoclonal antibodies and a fluorochrome labeled bacterial aerolysin (FLAER).<sup>2</sup> Through high sensitivity flow cytometry testing, as few as 0.01% PNH cells can be detected.<sup>1</sup>

## Clinical Indications

This test may be useful in the evaluation of patients with intravascular hemolysis, unexplained hemolysis, thrombosis with unusual features, or bone marrow failure.

## Results and/or Interpretation

Results are reported as:

- Negative
- Low-level PNH clone positive (.01 - <1%), with percentage
- PNH clone positive ( $\geq 1\%$ ), with percentage

## Limitations of the Assay

Results of PNH flow cytometry studies must be interpreted in the context of the clinical, laboratory and histologic findings.

Blood not collected in an EDTA tube or which has been improperly stored and handled prior to receipt cannot be processed.

Blood stability limit for PNH testing is 24 hours after stated draw time. Clinical significance of results on specimens 24-48 hours old should be evaluated in the context of other clinical and laboratory findings. Blood older than 48 hours from draw time cannot be analyzed.

## Methodology

High sensitivity flow cytometry for PNH is performed in accordance with Clinical Cytometry Society recommendations for high sensitivity flow cytometry testing.<sup>1</sup>

Cells are stained directly with FITC, PE, PE/Cy7, PerCP/Cy5.5, APC and APC/Cy7-labeled monoclonal antibodies to detect antigens of interest. Antibodies used include CD15, CD24, CD33, CD45, CD59 and glycophorin A. Additionally, FLAER is employed. For erythrocytes, antibodies to glycophorin A are used to specifically gate red cells and PNH clones are identified by lack of CD59 expression. PNH erythrocyte clones are divided into those with partial loss of CD59 (PNH Type II Cells) or complete loss of CD59 (PNH Type III Cells). For granulocytes, CD15 and CD33 are used to specifically gate granulocytes. PNH-type granulocytes are then identified by lack of expression of CD24 and lack of reactivity to FLAER.

## References

1. Borowitz MJ, Craig FE, DiGiuseppe JA *et al.* Guidelines for the Diagnosis and Monitoring of Paroxysmal Nocturnal Hemoglobinuria and Related Disorders by Flow Cytometry. *Cytometry B Clin Cytom.* 2010; 78B:211-230.
2. Sutherland DR, Kuek N, Davidson J *et al.* Diagnosing PNH with FLAER and Multiparameter Flow Cytometry. *Cytometry B Clin Cytom.* 2007; 72B: 167-177.

## Test Overview

<b>Test Name</b>	High Sensitivity Flow Cytometry for Paroxysmal Nocturnal Hemoglobinuria
<b>Reference Range</b>	Negative: No PNH clone detected
<b>Specimen Requirements</b>	Volume/Size: 4 mL; Type, Whole Blood; Tube/Container, EDTA (Lavender); Transport Temperature, Room temperature <b>Note:</b> Samples should be received in testing lab within 24 hours of stated draw time.
<b>Test Ordering Information</b>	Clearly indicate specimen source on sample label
<b>Billing Code</b>	81442
<b>CPT Codes</b>	88184, 88185, 88187, 86356

### Technical Information Contact:

Debbi Katanik, MT  
216.444.0042  
katanid@ccf.org

### Scientific Information Contact:

Eric D. Hsi, MD  
216.444.5230  
hsie@ccf.org